

NEW PRODUCT BULLETIN

Thermostrip Forehead Thermometer

Hyperthermia and Hypothermia Monitoring

- Transport through ER
- Easy to Read Moving Line
- Non-Invasive—Cost Effective
- Comparable to Esophageal Probes
- Core Adjusted Temperature Trends
- Fast, Disposable, Latex-Free

Continuous Display

A fast, simple temperature monitoring system based on proprietary color change technology.

To Read

1. Remove backing. Apply to center of forehead. Skin surface should be clean and dry.
2. Read temperature when yellow line stops moving. Typically 10-15 seconds.

Please see package insert for complete Instructions for Use.



PRODUCTS	
	01-0314
Description	Disposable (with adhesive backing) 50 thermometers per dispenser
°F or °C	°F & °C
Range	84°F-106°F / 29°C-41°C
Case	1

Liquid crystals are chemical compounds and mixtures that exhibit the mechanical properties of liquids and the optical characteristics of solids. Liquid crystals that change color with changes in temperature are known as thermochromic liquid crystals (TLC's). Scientifically they are identified as Cholesteric Liquid Crystals.

Thermostrip products utilize this exciting technology (TLC's) for the measurement of clinical temperature trends in the healthcare environment. Thermostrip products feature a continuous display of core-adjusted temperature as indicated by a bright line, which continuously moves across the temperature range.

www.thermometersite.com Phone: 800.585.9983 email: info@thermometersite.com

ThermometerSite.com is the premier location for temperature indicating color changing thermometers, thermographic supplies and reference material on the web.

We offer liquid crystal materials for non destructive research, testing and education.